



COVID-19 (March 2020) Infection Control Totes- Replenish as Required

Item	Description
1 Large Tote with Cover	To hold used items Prior to cleaning See step #1 in Table 1
1 Pair Heavy Rubber Gloves	See step #2, 11 in Table 1
Basin #1	See step #3 in Table 1
1 Bottle Pine-Sol Detergent	See step #4 in Table 1
1 Scrub Brush	See step #5 in Table 1
Basin #2	See step #8, 9, 10 in Table 1
2 Measuring Cups	See Table 3 & 4: Preparing Household Bleach as a Disinfectant & Dilution Rate Table
Ziploc Bags	Place dry, clean & disinfected equipment in a Ziploc bag until needed. Discard Ziploc bag after use. Do not reuse Ziploc bags



Infection Control Totes- Replenish as Required

Instruction Manual	Infection Prevention & Control Info and Procedures.
MSDS Sheets 3-Bleach-Sodium Hypochloride 3-Original Pin-Sol	Found in back of Blue Duo-Tang